

# App Guide

# Gait Phase Feedback 1

SageMotion  
Wearable Biofeedback System





## **Table of Contents**

Components	<i>(page 1)</i>
Wirelessly Connect to Computer or Cellphone	<i>(page 2)</i>
Gait Phase Feedback 1 App	<i>(page 3)</i>

# Components



Hub



Nodes (8x)



Battery



Node Straps: *Medium (8x), Short (4x), Long (2x)*




Cable A (10x)

-Connect Hub to Battery  
-Charge Nodes & Battery



Cable B (*optional use*)

-Connect Hub to Computer




Node Charging Station

# Wirelessly Connect to Computer or Cellphone

## 1) Connect Cable A to Battery and to Hub




## 2) On Computer/Cellphone, Connect to Wi-Fi: "SageMotion"



*Note 1: Need to wait for up to 1 minute for "SageMotion" to appear in Wi-Fi list. If it doesn't appear, try turning the Wi-Fi off and then on again on the computer/cellphone.*

*Note 2: Hub is connected after clicking "Connect" even if in Windows it shows "Connecting" or "No internet, open".*

## 3) On Computer/Cellphone, in Chrome Address Bar, Go To <http://192.168.12.1>



**[Note] If Computer Doesn't Have Wi-Fi:** *plug in Cable B to the Hub and to the ethernet port of your computer, then in chrome address bar, go to **<http://192.168.137.1>***

# Gait Phase Feedback 1 App

*The purpose of the Gait Phase Feedback 1 App is to measure gait phase of both feet and provide feedback during stance.*

## 1) Turn on 4 Nodes




*Slide switch toward middle to turn node on*



*Green light will blink after the node is on and running*

## 2) Select “Gait Phase Feedback 1” App



## 3) Click “Search”

Node List







# Gait Phase Feedback 1 App (*cont.*)

4) Configure 2 Sensor Nodes and 2 Feedback Nodes as Shown Below:

Node List

Search

Connect

Type	Position	MAC	
sensor	foot_left	88:6B:0F:E1:D8:9E	
sensor	foot_right	88:6B:0F:E1:D8:A6	
feedback	shank_left	88:6B:0F:E1:D8:9F	
feedback	shank_right	88:6B:0F:E1:D8:96	

5) Click “Connect”

Node List

Search

Connect

6) “Ready to collect data” Will Appear after Node Connection is Complete

Gait Phase Feedback 1 ▼

Start

Stop

✓ *Ready to collect data*



# Gait Phase Feedback 1 App (cont.)

7) Thread Straps through Nodes and Attach at Locations Shown Below:






*For both foot sensor nodes, the on/off switch points away from the body*

**[Note]** It is recommended to place the feedback nodes at the ankles as shown in the figure above, but they can be placed at any desired location

## How to Thread Straps



8) Click “Blink” for each Node to Confirm Correct Locations (red LED for given node blinks 3 times on click)

Type	Position	MAC			
sensor	foot_left	88:6B:0F:E1:D8:9E			
sensor	foot_right	88:6B:0F:E1:D8:A6			
feedback	shank_left	88:6B:0F:E1:D8:9F			
feedback	shank_right	88:6B:0F:E1:D8:96			



# Gait Phase Feedback 1 App (*cont.*)

## 9) In App Configuration, Enter Settings (Example Below)

### App Configuration

Trial Name	<input type="text" value="trial_1"/>
<b>Feedback Setting</b>	
Left Foot Feedback On?	<input type="text" value="true"/> ▼
Right Foot Feedback On?	<input type="text" value="true"/> ▼
When to Give Feedback?	<input type="text" value="Early stance"/> ▼
Feedback Pulse Length (sec)	<input type="text" value="0.5"/>
<b>Save Options</b>	
Save Mode	<input type="text" value="xlsx"/> ▼

# Gait Phase Feedback 1 App (*cont.*)

10) Click “Start” to Start Running the App



11) After the Trial is Finished, Click “Stop”



12) After Clicking “Stop”, a File from that Trial will Appear under Download Data. Click the File (e.g. trial1) to Download it to the Computer or Phone.

## Data Management

<input type="checkbox"/>	Name	Date▲	Duration	App	Type	Size	Rename	Delete
<input type="checkbox"/>	<u>trial1</u>	2021-12-12-11-39-18	0:00:18	Gait Phase Feedback 1	.xlsx	2.3 MB		

# Gait Phase Feedback 1 App (*cont.*)

## Description of Data in Downloaded File

**time** (sec): time since trial start

**Step\_Count\_Left**: left foot steps

**Gait\_Phase\_Left** : gait phase of left foot. 0 is “Early stance”; 1 is “Middle stance” ; 2 is “Late stance”; 3 is “Swing”

**Step\_Count\_Right**: right foot steps

**Gait\_Phase\_Right**: gait phase of right foot . 0 is “Early stance”; 1 is “Middle stance” ; 2 is “Late stance”; 3 is “Swing”

**Feedback\_Left**: feedback status for left foot. 0 is “feedback off”; 1 is “feedback on”

**Feedback\_Right** : feedback status for right foot. 0 is “feedback off”; 1 is “feedback on”

**SensorIndex\_1/2**: index of raw sensor data

**AccelX/Y/Z\_1/2** (m/s<sup>2</sup>): raw acceleration data

**GyroX/Y/Z\_1/2** (deg/s): raw gyroscope data

**MagX/Y/Z\_1/2** (μT): raw magnetometer data

**Quat1/2/3/4\_1/2** : quaternion data

**Sampletime\_1/2**: timestamp of each sensor

**Package\_1/2**: package number of each sensor